

Restriction Enzyme Troubleshooting Guide

The following guide can be used for troubleshooting restriction enzyme digestions. You may also be interested in reviewing additional tips for [optimizing digestion reactions](#).

Why is My Restriction Enzyme Not Cutting DNA?



Not getting the cleavage you expected? Let an NEB scientist help you troubleshoot your reaction.

| Problem | Possible Cause | Solution |
|-----------------------------|---|--|
| Incomplete or no digestion | Enzyme is inactive | <ul style="list-style-type: none"> ▪ Test enzyme on control DNA with known multiple sites ▪ Enzyme should be stored at -20°C. Enzymes stored at -70°C will freeze, and repeated thaw/freeze cycles may reduce enzyme activity |
| | Reaction conditions are not optimal | <ul style="list-style-type: none"> ▪ Use recommended buffer supplied with restriction enzyme ▪ Follow double digest recommendations, or try a sequential digest ▪ Repeat with fresh buffer. Additives present in buffer (e.g., DTT, SAM) may degrade over time |
| | Enzyme concentration is too low | <ul style="list-style-type: none"> ▪ Some plasmids or genomic DNAs may require up to 10–20 units/μg |
| | Additive is missing | <ul style="list-style-type: none"> ▪ Repeat reaction setup, being sure that enzyme and/or additives (e.g., BSA) are added |
| | DNA concentration is not optimal | <ul style="list-style-type: none"> ▪ NEB recommends 1 μg of DNA in a 50 μl reaction. ▪ Excess DNA may result in incomplete cleavage. |
| | Incubation time was too short | <ul style="list-style-type: none"> ▪ Some enzymes can exhibit slower cleavage towards specific sites. In most cases, 1–2 hours are sufficient. |
| | DNA is contaminated with an inhibitor | <ul style="list-style-type: none"> ▪ Assay substrate DNA in the presence of a control DNA. Control DNA will not cleave if there is an inhibitor present. Miniprep DNA is particularly susceptible to contaminants. ▪ Clean DNA with a spin column, resin or drop dialysis, or increase volume to dilute contaminant. |
| | Recognition site is not present | <ul style="list-style-type: none"> ▪ Confirm DNA sequence |
| | Cleavage is blocked by methylation | <ul style="list-style-type: none"> ▪ DNA isolated from a bacterial source may be blocked by Dam and Dcm methylation. DNA should be passed through a <i>dam-dcm</i>- strain (NEB #C2925). ▪ Eukaryotic genomic DNA may be blocked by CpG methylation. This can be overcome by cloning into a bacterial host. ▪ PCR products are not methylated |
| | DNA may be supercoiled | <ul style="list-style-type: none"> ▪ Restriction enzymes cleave supercoiled DNA with varying efficiency. Additional enzyme may be required. |
| Recognition site may be too | <ul style="list-style-type: none"> ▪ As a general rule, add 6 base pairs on either side of the recognition site for efficient cleavage | |

| | | |
|-----------------------------|--|---|
| | close to the end of the DNA fragment | |
| | Site preference | ▪ Enzyme requires two recognition sites for efficient cleavage. |
| | Quality of DNA | ▪ DNA purification procedures that use spin columns can result in high salt levels which inhibits enzyme activity. To prevent this, DNA solution should be no more than 25% of reaction volume. |
| Unexpected Cleavage Pattern | Determine nature of pattern | ▪ Run uncut substrate DNA alongside the digest. A partial digest will show bands found in the uncut, whereas star activity will show bands of unexpected size. |
| | DNA sample is contaminated | ▪ Prepare a new DNA sample |
| | Additional recognition sites are present in DNA | ▪ Confirm DNA sequence |
| | Star Activity | ▪ See tips for avoiding star activity (see page 5) and/or use a High-Fidelity (HF™) Restriction Enzyme (see page 4). |
| Smearing of DNA on gel | Enzyme has a high binding affinity to DNA and will not dissociate well | ▪ Add SDS to the gel loading dye/stop solution to a final concentration of 0.1–0.5% to help dissociate the enzyme from the DNA, or treat with protease before loading |
| | Nuclease contamination | ▪ Care should be taken to avoid cross contamination when setting up reactions |
| | Agarose running conditions | ▪ Use fresh running buffer and appropriate voltage to avoid overheating |
| | Master Mix containing dye | ▪ Use fresh gel |
| | Agarose gel or running buffer is old when using RE-Mix®. | ▪ Use fresh buffer, a fresh gel or add stop dye with EDTA |

HF™ is a trademark of New England Biolabs, Inc.

RE-Mix® is a registered trademark of New England Biolabs, Inc.